

FEATURE	RESTEP DVT System	Huntleigh* Flowtron Universal	CTC* Vasopress	Covidien* SCD Express
INDICATIONS				
Prevention of Deep Vein Thrombosis	Yes	Yes	Yes	Yes
CONSTRUCTION				
PVC Free	Yes	No	No	No
Thickness	0.876mm	4.69mm	8.81mm	0.927mm
Weight	58.9g	68.4g	107g	105.2g
PUMP FEATURES				
Battery "run" time	8 hours	4 hours	N/A	6 hours
Pump Weight	1.5lb	8.5lb	6.4lb	3.5lb
Power (max)	9 Watt	25 Watt	15 Watt	50 Watt
SAFETY/ALARMS				
High Pressure	Yes	Yes	No	No
Garment NOT deflating	Yes	Yes	No	No

*Information acquired from respective manufacturer websites.

Stryker Sustainability Solutions is the leading provider of reprocessing (cleaning, testing/verifying, packaging and sterilizing) and remanufacturing (disassembling, repairing and manufacturing) services for medical devices as well as comprehensive recycling and redistribution initiatives. Stryker is the only global OM (original manufacturer) to offer innovative, cost-effective reprocessing and remanufacturing programs that help meet the demands of its hospital partners. Specifically, Stryker Sustainability Solutions' programs target inefficient practices and reclaim and redirect consumable and financial resources in order to promote healthy and responsible hospitals. Stryker Sustainability Solutions' customers include most of the "Honor Roll" hospitals recognized by *U.S. News & World Report*. For more information about these programs, visit restep.stryker.com.